

AroCast

**PARAMOUNT**[®]

PARAMOUNT SURGIMED LTD.

AroCast

Made of pure Polyster Fiber and Synthetic Resin AroCast, Orthopedic tape, is better than existing Orthopedic plaster bandage for casting. AroCast effectively hardens and fixes the surface of the cast, and it works especially well for Orthopedic casting, cast bracing, weight bracing, replacement cast, etc.

STORAGE

AroCast should be kept at room temperature. Do not use the content if the package is damaged. Prolonged storage at high temperature and / or under the direct sunlight may reduce its quality and life-cycle.

AroCast is

Since AroCast is Orthopedic tape made of Polyester Fibers and Resins, it has the following properties.

Feature

- ✓ Dries fast
- ✓ Easy and simple to use.
- ✓ Light weight
- ✓ Vents Superbly
- ✓ Impact resistant.
- ✓ Easily applicable to any body area.
- ✓ Provides excellent flexibility and molding.
- ✓ Environmentally safe (burns up completely)



PARAMOUNT SURGIMED LTD.

1, L.S.C., Okhla Indl. Area, Phase-II, New Delhi -29 (India)
Tel : +91-11-26389812, 26389813, 41616333
Fax : +91-11-26389815, 41616555, 41621000
E-mail : psk@paramountblades.com
Website : <http://www.paramountblades.com>