



**PARAMOUNT**  
Total Blades Solution Company

# Disposable Scalpels

## disposable scalpels

Consist of pre - assembled stainless steel surgical blades mounted on ABS plastic handle with fluted finger grips and the blades are safely protected by LDPE guard. These are individually packed in peel-open paper/poly pouches and 10 such pcs are packed in box of 10's. Color of handles can be changed to customer's requirements.



## stitch cutters

The long and short stitch cutters require no handles. The long version is made of stainless steel and the short version is of carbon steel. The mini stainless steel stitch cutter can be used with handle Nos. 3L, 7 & 9. Alternatively, it is also available on a pre-mounted ABS disposable plastic handle.



**PARAMOUNT SURGIMED LTD.**

1, L.S.C., Okhla Indl. Area, Phase-II, New Delhi -20 (India)  
Tel. : +91-11-26389812, 26389813, 41616333  
Fax : +91-11-26389815, 41616555, 41021000  
E-mail : [psl@paramountblades.com](mailto:psl@paramountblades.com)  
Website : <http://www.paramountblades.com>