


PARAMOUNT®

Total Blades Solution Company

IS : 3319




CM/L 8536384



surgical blades

Available in Stainless Steel and Carbon Steel. Each blade is protected by VCI paper in a peel-open foil sachet. It is Sterilized by Gamma Radiation of 2.5 M Rads.




surgical skin graft knife blades

Surgical Skin Graft Blades are 158 mm long, individually wrapped in VCI paper in single peel packets in packs of 10. These blades can be used with Watson, Braithwaite or Cobbett Knives.

blood lancet

supplied sterile and individually packed. lancet comes in a pack of 200 pieces per box. the tip is specially designed for easy pricking, which minimizes patient trauma.



metal handles & postmortem blades

Metal Handles made of Stainless Steel, autoclavable and are available in sizes 3, 4, 3L, 4L, 5, 7, 8, 9. Postmortem blades are available in sizes 40, 40B, 60, 60B, 61, & 62.

